

OctolineTM₂

8-in-1 one-stop multi skin care system. Makes healthy and shiny skin by applying ultrasound, Iontophoresis, and air jet-peel function.

Multi Programs

- Various Program for skin type: Aging skin, acne skin, Sensitive skin
- 4-modes: UltraSon mode, Sono mode, Kion mode and jet mode

STRENGTHS

- Various Skin Care
- Simultaneous use of two functions
- Efficiency, stability and economic feasibility
- Wide operation panel for user convenience
- Simple and convenient operation




OctolineTM₂

Jet Peel & Spray effectively exfoliates the skin with powerful air pressure and purified water without damaging the epidermis. It improves skin tone and skin texture.

JET MODE

- Dead Skin cell peeling and skin regeneration effect
- Jet tip and Spray




Jet / Spray

OctolineTM₂

Ultrasonic scrubber removes dead skin cells remaining on the epidermis through low stimulation micro vibrations of 20,000 to 30,000 times per second and also removes sebum and wastes in pores.

SONO MODE

- Deep Cleansing effect
- Cleansing, Scaling, Lifting, Vitamin
- Adjustable 5 levels






Scrubber

Iontophoresis is a method of rapidly penetrating ionized solution using electric current. It quickly and safely penetrates active ingredients and allows optimal treatment in various areas using various types of handpieces.

ION MODE

- Effective solution infiltration and muscles massage effect
- +Ion, -Ion, Vitamin, Lifting
- Pin, roll Y-type handpiece
- Adjustable 10 Levels



Iontophoresis

OctolineTM₂

Heat generated by ultrasound relaxes the skin layer and penetrates the solution quickly and effectively. In particular, solutions with large molecular sizes can penetrate, which is effective in nourishing.



ULTRASON MODE

- Massage effect by the vibration
- Wrinkle, Acne, Couperose, Cellulite
- Handpieces for face and dark circle
- Adjustable 5 levels



Ultrasound